

### **Ethanol precipitation of DNA.**

Measure plasmid volume. Bring to 200  $\mu$ l with dH<sub>2</sub>O.

Add 20  $\mu$ l of 3M sodium acetate (pH 5.2). Mix by vortexing.

Add 660  $\mu$ l of cold 100% EtOH (keep in your -20°C).

Mix and put at -80°C for at least 30 min. Good stopping point.

Spin at max speed in microfuge in cold room for 15 min.

Carefully decant supernatant. Add 1 ml of 70% EtOH at room temp. Spin 3 min at full speed in microfuge. Decant carefully. Use pipette to remove remaining ethanol.

Briefly airdry and resuspend in 30  $\mu$ l of nuclease-free water. Nanodrop.