

General protein expression protocol.

- 1) In morning (around 9-10am), pick single colony of expression strain into 5 ml of LB+appropriate antibiotics and shake at 37°C.
- 2) At end of day (around 6 pm) transfer 5 ml culture into 1L of LB+appropriate antibiotics and shake overnight at 37°C.
- 3) The next morning check the OD600. You are aiming for 0.8-1.0, unless otherwise indicated. When culture reaches this OD600, take an **uninduced sample**.
 - take 500 µl of culture and pellet in a microfuge (fullspeed, 30-60 sec)
 - aspirate of supernatant
 - resuspend in 100 µl of 1x protein sample buffer
 - boil 5 minutes
 - keep on ice
- 4) Induce culture with indicated amount of IPTG. Shake for 4 hours at the desired temperature. I use 30 °C for MBP tagged proteins and 16 °C for most other preps.
- 5) **Take induced sample.**
 - take 500 µl of culture and pellet in a microfuge (fullspeed, 30-60 sec)
 - aspirate of supernatant
 - resuspend in 200 µl of 1x protein sample buffer
 - boil 5 minutes
 - keep on ice
- 6) **Take soluble sample.**
 - take 10 ml of culture and pellet in centrifuge (4000 rpm, 10 min, 4°C)
 - aspirate of supernatant
 - freeze until ready to process (see step 8)
- 7) Pellet 1 L culture in centrifuge (4000 rpm, 15 min, 4°C). Wash with 25 mls of PBS and resuspend culture. Transfer to a 50 ml conical tube and bring volume up to 50 ml with PBS. Invert several times to mix the sample. Pellet in centrifuge (4000 rpm, 10 min, 4°C). Freeze in liquid nitrogen and store at -20°C.
- 8) Processing soluble samples. Resuspend in 1ml PBS+Complete protease inhibitor (use 1 Mini tablets per 5 ml) and transfer to a 1.5 ml Eppendorf, keeping sample on ice.
 - lyse the bacteria by sonication keeping all reagents on ice during sonication: (Sonication: 2x30 sec, 15% amplitude. 1x30 sec 25% amplitude)
 - pellet insoluble debris by centrifugation (14k, 15 min, 4°C).
 - make soluble input fractions by taking 100 µl of lysate and adding 30 µl of 4xSample
 - store samples at -20°C